

General suggestion

Always do this as a team. One person performs the alignment and one person monitor the actions of the first person.

- 1.) Turn off X-rays
- 2.) Move detector back
- 3.) Move LT out of the way
- 4.) Remove beam stop
- 5.) Place PIN diode on collimator and secure with black tape
- 6.) Plug PIN diode powder supply into the outlet
- 7.) Plug the dual plug into multi-meter
- 8.) Turn dial on multi-meter to 200mv DC volts (meter will read 0.0)
- 9.) Place Radiation meter below collimator on platform and turn on
- 10.) With the service key override the safety interlocks
- 11.) Turn X-rays on move X-ray generator Key to 2nd position
- 12.) Bring X-rays to 40 kV and 40mA
- 13.) On the SMART or APEX instrument: start SMART (or APEX)
- 14.) Shut the X-ray enclosure doors
- 15.) Open safety shutter
 - a. On D8 push and hold red shutter button on shutter control box (just below the doors)
 - b. On SMART1000 push the shutter open button on shutter control box (this is the reset button in normal mode)
- 16.) Safety shutter will now open
- 17.) Check for stray X-rays
 - a. GM meter should be silent (except for normal background radiation)
- 18.) On the SMART or APEX: toggle rotary shutter with SMART (or APEX)
- 19.) The GM meter should read a positive MV reading
 - a. (negative if dual plug is plugged in wrong)

For a new X-ray tube and a 0.5mm collimator (only front pinhole in)

Mo Ceramic Tube	40kV/40ma	143 mV
Mo Glass Tube	40kV/40ma	120 mV
Cu Ceramic Tube	40kV/40ma	90 mV
Cu Glass Tube	40kV/40ma	50 mV
- 20.) Check for stray X-rays
 - a. If no stray radiation is seen then it is safe to open the enclosure doors
 - i. Only one person should have their hands inside the enclosure box!
 - ii. Only put your hands and arms into the enclosure box!
 - iii. Work quickly but safely
- 21.) Use a second GM to measure around the tube and monochromator
 - a. Do a through sweep of the area to measure for stray radiation
 - b. Never assume that it is safe!
- 22.) If no stray radiation is detected then you may begin.