

Practicals : INS preparation.

J. Reibenspies

V1.0\_2009

Data Preparation | XPREP

Command Line

```
Directory of c:\Structures\yus1
09/10/2008  03:10 PM                163,685 data0m.p4p
09/10/2008  03:10 PM                163,685 yus1.p4p
                2 File(s)            327,370 bytes
                0 Dir(s)  151,593,357,312 bytes free

c:\Structures\yus1>dir *.hkl
Volume in drive C is OS
Volume Serial Number is 6EDB-0B0A

Directory of c:\Structures\yus1
09/10/2008  04:12 PM                821,760 yus1.hkl
                1 File(s)            821,760 bytes
                0 Dir(s)  151,593,357,312 bytes free

c:\Structures\yus1>
c:\Structures\yus1>
c:\Structures\yus1>
c:\Structures\yus1>xprep yus1 data0m.p4p
```

Signal / Noise

```
27391 Reflections read from file yus1.hkl
Mean (I/sigma) = 8.56

Lattice exceptions: P      A      B      C      I      F      Obv      Rev      All
N (total) =           0 13672 13679 13671 13680 20511 18269 18307 27391
N (int>3sigma) =       0 10090 10093 10127  9955 15155 13422 13487 20253
Mean intensity =    0.0 28.5 25.8 29.0 28.6 27.7 29.0 29.0 28.9
Mean int/sigma =    0.0  8.5  8.5  8.6  8.5  8.5  8.5  8.6  8.6

Select option [P]: █
```

Select option [H]:

Determination of reduced (Niggli) cell

Transformation from original cell (HKLf-matrix):

```
-1.0000  0.0000  0.0000  0.0000 -1.0000  0.0000  0.0000  0.0000  1.0000
```

Unitcell: 8.661 16.319 16.425 89.97 89.96 90.00

Niggli form: a.a = 75.01 b.b = 266.32 c.c = 269.77  
 b.c = 0.12 a.c = 0.11 a.b = 0.01

Search for higher METRIC symmetry

Identical indices and Friedel opposites combined before calculating R(sym)

---

Option A: FOM = 0.184 deg. TETRAGONAL P-lattice R(sym) = 0.544 [ 4685]  
Cell: 16.319 16.425 8.661 89.96 90.00 89.97 Volume: 2321.41  
Matrix: 0.0000 -1.0000 0.0000 0.0000 0.0000 1.0000 -1.0000 0.0000 0.0000

Option B: FOM = 0.051 deg. ORTHORHOMBIC P-lattice R(sym) = 0.077 [ 4537]  
Cell: 8.661 16.319 16.425 90.03 90.04 90.00 Volume: 2321.41  
Matrix: 1.0000 0.0000 0.0000 0.0000 1.0000 0.0000 0.0000 0.0000 1.0000

---

Option C: FOM = 0.025 deg. MONOCLINIC P-lattice R(sym) = 0.060 [ 2891]  
Cell: 8.661 16.319 16.425 90.03 90.04 90.00 Volume: 2321.41  
Matrix: 1.0000 0.0000 0.0000 0.0000 1.0000 0.0000 0.0000 0.0000 1.0000

---

Option D: FOM = 0.044 deg. MONOCLINIC P-lattice R(sym) = 0.062 [ 2795]  
Cell: 16.319 8.661 16.425 90.04 90.03 90.00 Volume: 2321.41  
Matrix: 0.0000 -1.0000 0.0000 -1.0000 0.0000 0.0000 0.0000 0.0000 -1.0000

---

Option E: FOM = 0.051 deg. MONOCLINIC P-lattice R(sym) = 0.065 [ 2892]  
Cell: 8.661 16.425 16.319 89.97 90.00 90.04 Volume: 2321.41  
Matrix: 1.0000 0.0000 0.0000 0.0000 0.0000 1.0000 0.0000 -1.0000 0.0000

---

Option F retains original cell

Select option [B]:

```

Current dataset: yus1.hkl                      Wavelength: 1.54178 Chiral: ?
-----
Original cell:  8.661  16.319  16.425  90.03  90.04  90.00  Vol  2321.4
                Esds:  0.003  0.005  0.005  0.03  0.03  0.02  Lattice: P
-----
Current cell:   8.661  16.319  16.425  90.03  90.04  90.00  Vol  2321.4
-----
Matrix:  1.0000  0.0000  0.0000  0.0000  1.0000  0.0000  0.0000  0.0000  1.0000
-----
Crystal system: Orthorhombic                  Lattice: P
-----

[D] Read, modify or merge DATASETS           [C] Define unit-cell CONTENTS
[P] Contour PATTERSON sections                [F] Set up shelxtl FILES
[H] Search for HIGHER metric symmetry        [R] RECIPROCAL space displays
[S] Determine or input SPACE GROUP           [U] UNIT-CELL transformations
[A] Absorption, powder, SIR, SAD, MAD etc.   [T] Change TOLERANCES
[M] Test for MEROHEDRAL TWINNING            [O] Self-rotation function
[L] Reset LATTICE type of original cell     [Q] QUIT program

Select option [S]: █

```

```

[S] Determine SPACE GROUP
[C] Must be CHIRAL (sample is optically active)
[N] NOT NECESSARILY chiral (eg. may be racemate)
[I] INPUT known space group
[E] EXIT to main menu or [Q] QUIT program

Select option [S]:

[A] Triclinic, [M] Monoclinic, [O] Orthorhombic, [T] Tetragonal,
[H] Trigonal/Hexagonal, [C] Cubic or [E] EXIT

Select option [O]:

Lattice exceptions: P      A      B      C      I      F      Obv      Rev      All
N (total) =              0 13672 13679 13671 13680 20511 18269 18307 27391
N (int>3sigma) =         0 10090 10093 10127  9955 15155 13422 13487 20253
Mean intensity =         0.0 28.5 25.8 29.0 28.6 27.7 29.0 29.0 28.9
Mean int/sigma =         0.0  8.5  8.5  8.6  8.5  8.5  8.5  8.6  8.6

Lattice type [P, A, B, C, I, F, O(obv.), R(rev. rhomb. on hex. axes)]

Select option [P]: █

```

Select option [S]:

[A] Triclinic, [M] Monoclinic, [O] Orthorhombic, [T] Tetragonal,  
[H] Trigonal/Hexagonal, [C] Cubic or [E] EXIT

Select option [O]:

Lattice exceptions:	P	A	B	C	I	F	Obv	Rev	All
N (total) =	0	13672	13679	13671	13680	20511	18269	18307	27391
N (int>3sigma) =	0	10090	10093	10127	9955	15155	13422	13487	20253
Mean intensity =	0.0	28.5	25.8	29.0	28.6	27.7	29.0	29.0	28.9
Mean int/sigma =	0.0	8.5	8.5	8.6	8.5	8.5	8.5	8.6	8.6

Lattice type [P, A, B, C, I, F, O(obv.), R(rev. rhomb. on hex. axes)]

Select option [P]:

Mean  $|E^*E-1| = 0.777$  [expected .968 centrosym and .736 non-centrosym]

Systematic absence exceptions:

	b--	c--	n--	21--	-c-	-a-	-n-	-21-	--a	--b	--n	--21
N	1315	1303	1308	35	787	806	783	62	753	749	752	94
N I>3s	846	863	777	0	496	542	494	4	566	534	478	16
<I>	33.9	33.8	24.1	0.3	35.4	27.8	29.1	0.2	40.7	73.0	60.3	0.5
<I/s>	8.0	7.8	7.0	0.5	8.2	8.4	7.5	0.9	9.5	8.8	7.5	1.2

Identical indices and Friedel opposites combined before calculating R(sym)

Option	Space Group	No.	Type	Axes	CSD	R(sym)	N(eq)	Syst. Abs.	CFOM
[A]	P2(1)2(1)2(1)	# 19	chiral	1	5917	0.077	4537	1.2 / 7.0	1.89

Select option [A]:

```
[D] Read, modify or merge DATASETS          [C] Define unit-cell CONTENTS
[P] Contour PATTERSON sections              [F] Set up shelxtl FILES
[H] Search for HIGHER metric symmetry      [R] RECIPROCAL space displays
[S] Determine or input SPACE GROUP         [U] UNIT-CELL transformations
[A] Absorption, powder, SIR, SAD, MAD etc. [T] Change TOLERANCES
[M] Test for MEROHEDRAL TWINNING          [O] Self-rotation function
[L] Reset LATTICE type of original cell    [Q] QUIT program
```

Select option [D]: █

```
[M] Sort-MERGE current data (no scaling)    [C] Change CURRENT dataset
[I] LEAST-SQUARES scale and merge datasets  [W] WRITE dataset to file
[I] INCLUDE Rfree flags from another file   [R] READ in another dataset
[S] Display intensity STATISTICS           [D] DELETE stored dataset
[F] FACE-indexed absorption corrections     [P] PSI-scan absorption corr.
[C] Copy file, TRANSFORM hkl and cosines   [A] MAD, SAD, SIR or SIRAS
[H] Apply HIGH/low resolution cutoffs      [N] NORMALIZE/scale sigmas
[G] Generate simulated powder diagrams     [U] Anisotropic scaling
[Z] Expand data to triclinic               [E] EXIT to main menu
[X] Parsons' Q values and Flack x parameter [Q] QUIT program
```

```
[R] Define resolution ranges (currently selected automatically)
[P] Output R(pim) rather than R(sigma)
[N] Do NOT merge reflections
[I] Merge IDENTICAL indices only
[S] Merge SYMMETRY equivalents (but not Friedel opposites)
[A] Merge ALL equivalents (including Friedel opposites)
[E] EXIT to merge datasets menu
[Q] QUIT program
```

Resolution	#Data	#Theory	%Complete	Redundancy	Mean I	Mean I/s	Rint	Rsigma
Inf - 2.60	106	106	100.0	28.51	100.8	65.82	0.0546	0.0156
2.60 - 2.00	111	111	100.0	28.93	57.7	66.39	0.0533	0.0138
2.00 - 1.70	113	113	100.0	24.50	35.9	53.40	0.0529	0.0151
1.70 - 1.54	107	107	100.0	25.19	21.8	48.21	0.0826	0.0179
1.54 - 1.42	113	113	100.0	25.57	21.3	47.12	0.0941	0.0177
1.42 - 1.33	110	110	100.0	19.11	16.4	36.08	0.1585	0.0250
1.33 - 1.25	121	121	100.0	13.80	15.3	28.47	0.1522	0.0292
1.25 - 1.19	119	119	100.0	11.47	15.2	25.47	0.1268	0.0387
1.19 - 1.14	123	123	100.0	11.03	13.3	23.76	0.1620	0.0398
1.14 - 1.09	135	135	100.0	10.55	9.2	19.01	0.1848	0.0544
1.09 - 1.05	131	131	100.0	9.40	7.4	13.23	0.2651	0.0690
1.05 - 1.02	116	116	100.0	5.53	4.8	6.72	0.3508	0.1539
1.02 - 0.99	123	123	100.0	4.95	5.5	7.21	0.2937	0.1513
0.99 - 0.96	140	140	100.0	4.55	4.4	6.06	0.3166	0.2024
0.96 - 0.93	157	157	100.0	4.55	4.1	5.88	0.3621	0.2036
0.93 - 0.91	129	129	100.0	4.19	4.5	6.36	0.3140	0.1626
0.91 - 0.89	115	115	100.0	3.97	2.9	4.60	0.3884	0.2292
0.89 - 0.87	28	39	71.8	1.44	2.0	2.63	0.4537	0.7050
-----								
0.97 - 0.87	532	543	98.0	4.09	3.8	5.54	0.3461	0.2167
Inf - 0.87	2097	2108	99.5	12.99	18.4	25.46	0.0790	0.0369

Merged [A], lowest resolution = 16.42 Angstroms, 1709 outliers downweighted

Enter <CR> to continue

- |  |                               |
|--|-------------------------------|
| [D] Read, modify or merge DATASETS         | [C] Define unit-cell CONTENTS |
| [P] Contour PATTERSON sections             | [F] Set up shelxtl FILES      |
| [H] Search for HIGHER metric symmetry      | [R] RECIPROCAL space displays |
| [S] Determine or input SPACE GROUP         | [U] UNIT-CELL transformations |
| [A] Absorption, powder, SIR, SAD, MAD etc. | [T] Change TOLERANCES         |
| [M] Test for MEROHEDRAL TWINNING           | [O] Self-rotation function    |
| [L] Reset LATTICE type of original cell    | [Q] QUIT program              |

Select option [C]:

Select option [E]: F

Enter formula; numbers follow elements or brackets, 2nd character of element name must be lower case, may include: Me, Et, Ep, Bu, Ph, or Cp:  
C20 H40 N4 O3 Cl2

Tentative Z (number of formula units/cell) = 4.0 giving rho = 1.303,  
non-H atomic volume = 20.0 and following cell contents and analysis:

C	80.00	52.74 %	H	160.00	8.85 %
N	16.00	12.30 %	O	12.00	10.54 %
Cl	8.00	15.57 %			

[Z] change Z, [F] new FORMULA, [R] change RADIATION,  
[E] EXIT to main menu or [Q] QUIT program

Select option [E]: █

Select option [F]:

Output file name (without extension) [yus1]:

XM/SHELXD (M) or XS/SHELXS (S) format [S]:

File yus1.ins set up as follows:

```
TITL yus1 in P2(1)2(1)2(1)
CELL 1.54178 8.6608 16.3192 16.4246 90.000 90.000 90.000
ZERR 4.00 0.0029 0.0046 0.0049 0.000 0.000 0.000
LATT -1
SYMM 0.5-X, -Y, 0.5+Z
SYMM -X, 0.5+Y, 0.5-Z
SYMM 0.5+X, 0.5-Y, -Z
SFAC C H N O CL
UNIT 80 160 16 12 8
TEMP 0
SIZE 0.01 0.07 0.13
TREF
HKLF 4
END
```

Do you wish to (over)write the intensity data file yus1.hkl ? [Y]: █

Output : \*.hkl \*.ins (INS for XS or XM)