


Practical Guide to Structure Solution, Refinement and Publication.

SHELXTL

Instructor : J. Reibenspies, Ph. D.
Nattamai Bhuvanesh, Ph.D.

Version 1.0.0

Flow Chart



Shell Programs

SHELXTL – BRUKER

- a) license stipulates that program can be loaded only on computers owned by the university.

XSEED – L. Barbour

- a) license stipulates that program must be used by university employees. If you leave the university then your license is no longer valid.

WINGX – L. Farrugia

- a) Freeware, must register with Farrugia.

Other programs

SHELXS, SHELXL - Free must register with G. Sheldrick

SIR88, SIR97, SIR2002 – Free must sign license with C. Giacovazzo

RASMOL, POVRAY, MERCURY

Data Input: 3I4, 2F8.2 $h, k, l, I, \sigma(I)$ [* .hkl]

Parameter file from SMART [* .p4p]

Information required :

- a) Molecular Formula (unit cell contents)
- b) Anticipated structure

Output crystallographic information file CIF

Crystal Data

Coordinates

Bond lengths and angles


Thermal parameters

Hydrogen atoms

Fo/Fc information FCF

Getting Started. SHELXTL

- Move the *.hkl and *.p4p file to a new directory
- Start SHELXTL
 - Point to PROJECT and then NEW.
 - Find your files and name the project to begin



- Point to XPREP and start program
 - A new window will appear

```

XPREP Ver 6.12 W95/98/NT/2000/ME Copyright Bruker-AXS 2001
+ XPREP - Reciprocal space exploration - Version 6.12 - W95/98/NT/2000/ME +
+ COPYRIGHT(c) 2001 Bruker-AXS All Rights Reserved +
+-----+
Screen size: 1024 x 768
Window size: 640 x 693
Font size: 8 x 16
Number of colors: 256

8226 Reflections read from file sucrose.hkl
Mean (I/sigma) = 5.55

Lattice exceptions: P A B C I F Obv Rev All
N (total) = 0 4117 4094 4095 4108 6153 5483 5485 8226
N (int>3sigma) = 0 2153 2168 2207 2175 3264 2933 2923 4397
Mean intensity = 0.0 18.2 18.2 16.3 17.5 17.6 16.9 17.1 17.0
Mean int/sigma = 0.0 5.5 5.7 5.7 5.7 5.7 5.6 5.7 5.6

Lattice type [P, A, B, C, I, F, O(obv.), R(rev. rhomb. on hex. axes)]
Select option [P]: █

```

in this case choose the primitive cell (based on the statistics)
Mean Intensity >> mean I/sigma(I) for A,B,C,I, F etc.
Number of I > 3 sigma(I) is large for A,B,C,I,F, etc.

Search for Higher Symmetry.

```

Select option [H]: H

Determination of reduced (Niggli) cell
Transformation from original cell (HKL-matrix):
-1.0000 0.0000 0.0000 0.0000 1.0000 0.0000 0.0000 0.0000 -1.0000
Unitcell: 7.758 8.702 10.861 90.00 102.95 90.01
Niggli form: a.a = 60.18 b.b = 75.73 c.c = 117.96
             b.c = -0.01 a.c = -18.88 a.b = -0.01

Search for higher METRIC symmetry
Identical indices and Friedel opposites combined before calculating R(sym)

-----
Option A: FOM = 0.013 deg. MONOCLINIC P-lattice R(sym) = 0.040 [ 1464]
Cell: 7.758 8.702 10.861 90.00 102.95 89.99 Volume: 714.58
Matrix: 1.0000 0.0000 0.0000 0.0000 1.0000 0.0000 0.0000 0.0000 1.0000
-----
Option B retains original cell

Select option [A]: █

```

In this case the Monoclinic Cell is chosen (notice an Rsym > .1 may indicate that the choice of the cell is wrong)

Select Determine or Select Space Group from main menu.

```

XPREP Ver 6.12 W95/98/NT/2000/ME Copyright Bruker-AXS 2001
Current dataset: sucrose.hkl          Wavelength: 0.71073 Chiral: ?
-----
Original cell:  7.758  8.702 10.861  90.00 102.95  89.99  Vol  714.6
                Esds:  0.002  0.002  0.002  0.00  0.00  0.00  Lattice: P
-----
Current cell:  7.758  8.702 10.861  90.00 102.95  89.99  Vol  714.6
-----
Matrix: 1.0000  0.0000  0.0000  0.0000  1.0000  0.0000  0.0000  0.0000  1.0000
-----
Crystal system: Monoclinic          Lattice: P
-----
[D] Read, modify or merge DATASETS          [C] Define unit-cell CONTENTS
[P] Contour PATTERSON sections              [F] Set up shelxtl FILES
[H] Search for HIGHER metric symmetry      [R] RECIPROCAL space displays
[S] Determine or input SPACE GROUP         [U] UNIT-CELL transformations
[A] Absorption, powder, SIR, SAD, MAD etc. [T] Change TOLERANCES
[M] Test for MEROHEDRAL TWINNING          [O] Self-rotation function
[L] Reset LATTICE type of original cell    [Q] QUIT program

Select option [S]: S

```

Now set up the automatic space group determination

```

Crystal system: Monoclinic          Lattice: P
-----
[S] Determine SPACE GROUP
[C] Must be CHIRAL (sample is optically active)
[N] NOT NECESSARILY chiral (eg. may be racemate)
[I] INPUT known space group
[E] EXIT to main menu or [Q] QUIT program

Select option [S]: s

[A] Triclinic, [M] Monoclinic, [O] Orthorhombic, [T] Tetragonal,
[H] Trigonal/Hexagonal, [C] Cubic or [E] EXIT

Select option [M]: M

Lattice exceptions: P      A      B      C      I      F      Obv      Rev      All
N (total) =              0  4117  4094  4095  4108  6153  5483  5485  8226
N (int>3sigma) =         0  2153  2168  2207  2175  3264  2933  2923  4397
Mean intensity =         0.0 18.2  18.2  16.3  17.5  17.6  16.9  17.1  17.0
Mean int/sigma =         0.0  5.5  5.7  5.7  5.7  5.7  5.6  5.7  5.6

Lattice type [P, A, B, C, I, F, O(obv.), R(rev. rhomb. on hex. axes)]

Select option [P]: P

```

```

Lattice exceptions: P      A      B      C      I      F      Obv      Rev      A1
N (total) =           0    4117    4094    4095    4108    6153    5483    5485    822
N (int>3sigma) =       0    2153    2168    2207    2175    3264    2933    2923    439
Mean intensity =     0.0    18.2    18.2    16.3    17.5    17.6    16.9    17.1    17.
Mean int/sigma =     0.0     5.5     5.7     5.7     5.7     5.7     5.6     5.7     5.

Lattice type [P, A, B, C, I, F, O(obv.), R(rev. rhomb. on hex. axes)]
Select option [P]: P

Mean |E*E-1| = 0.792 [expected .968 centrosym and .736 non-centrosym]

Systematic absence exceptions:

      -21-  -a-  -c-  -n-
N       18   304   301   297
N I>3s    0   164   156   154
<I>      0.1  37.6  35.4  36.9
<I/s>    0.4   6.8   6.0   6.5

Identical indices and Friedel opposites combined before calculating R(sym)

Option  Space Group  No.  Type  Axes  CSD  R(sym)  N(eq)  Syst. Abs.  CFOM
[A] P2(1)           # 4  chiral  1  3543  0.040  1464  0.4 / 5.6  2.11
[B] P2(1)/m        # 11 centro  1  402  0.040  1464  0.4 / 5.6  5.13

Select option [A]: █

```

notes :

The E^2-I statistics predict a non-centrosymmetric space group

$I > I/\sigma(I)$ for the a,c and n glides

$I < I/\sigma(I)$ for 2_1 screw axis

$N(I > 3 \sigma(I))$ is 0 for 2_1

Space group can be $P 2_1$ and $P 2_1/m$

$P 2_1/m$ is a centrosymmetric space group and not consistent with the E^2-I statistics.

The structure has a known chiral carbon center (i.e. the space group cannot contain a mirror plane or an inversion center)

Space $P 2_1$ is chosen.

(see the attached handout for help in choosing space groups)

```

[D] Read, modify or merge DATASETS           [C] Define unit-cell CONTENTS
[P] Contour PATTERSON sections               [F] Set up shelxtl FILES
[H] Search for HIGHER metric symmetry       [R] RECIPROCAL space displays
[S] Determine or input SPACE GROUP          [U] UNIT-CELL transformations
[A] Absorption, powder, SIR, SAD, MAD etc.  [T] Change TOLERANCES
[M] Test for MEROHEDRAL TWINNING           [O] Self-rotation function
[L] Reset LATTICE type of original cell     [Q] QUIT program

Select option [D]: C

```

Skip D (for now) and define unit-cell contents

```

Enter formula; numbers follow elements or brackets, 2nd character of element
name must be lower case, may include: Me, Et, Pp, Bu, Ph, or Cp:
C12 H22 O11

Tentative Z (number of formula units/cell) = 2.0 giving rho = 1.591,
non-H atomic volume = 15.5 and following cell contents and analysis:

C      24.00   42.10 %           H      44.00   6.48 %
O      22.00   51.42 %

[Z] change Z, [F] new FORMULA, [R] change RADIATION,
[E] EXIT to main menu or [Q] QUIT program

Select option [E]:

```

Enter the formula $C_{12}H_{22}O_{11}$ for sugar. Check the Tentative Z value (should be near 2 for this space group) and the density (see density handout for more information)
If ok choose E and return to main menu.

Choose F from the main menu and write the SHELXTL files.

```
Select option [D]: F

Output file name (without extension) [sucrose]:

File sucrose.ins set up as follows:

TITL sucrose in P2(1)
CELL 0.71073  7.7576  8.7022  10.8611  90.000  102.946  90.000
ZERR  2.00  0.0018  0.0020  0.0025  0.000  0.004  0.000
LATT -1
SYMM -X, 0.5+Y, -Z
SFAC C H O
UNIT 24 44 22
TEMP 0
TREF
HKLF 4
END

Do you wish to (over)write the intensity data file sucrose.hkl ? [N]: █
```

You may need to write a new file/project name if the program has transformed the unit cell. If you have a backup copy of the HKL file then overwrite the old file here.

Quit the program now and go back to the SHELXTL program

Point to XS to get started (If you had to generate a new file in xprep then go back to project and start a new project)

The XS window will start

```
C:\WINDOWS\System32\cmd.exe
+ XS - CRYSTAL STRUCTURE SOLUTION - SHELXTL Ver. 6.12 W95/98/NT/2000/ME +
+ Copyright(c) 2001 Bruker AXS All Rights Reserved +
+ sucrose started at 09:27:09 on 06-Jul-2003 +
+*****+
Read instructions and process reflection data
Data: 1803 unique, 1298 observed R(int) = 0.0571 R(sigma) = 0.0586
Systematic absence violations: 0 Bad equivalents: 0
ESEL Emin 1.200 Emax 5.000 Delu 0.005 renorm 0.700 axis 0
OMIT s 4.00 2theta(lim) 180.0
INIT nn 12 nf 16 s+ 0.800 s- 0.200 wr 0.200
PHAN steps 10 cool 0.900 Boltz 0.300 ns 181 mtpv 40 mnqr 10
TREF np 256. nE 252 kpscals 0.800 ntan 3 wn -0.750
FMAP code 8
PLAN npeaks -30 del1 0.500 del2 1.500
MORE verbosity 1
TIME t 9999999.

181 Reflections and 2190. unique TPR for phase annealing
227 Phases refined using 3690. unique TPR
227 Reflections and 3690. unique TPR for R(alpha)
1223 Unique negative quartets found, 1223 used for phase refinement
352 Unique NQR employed in phase annealing
128 Parallel refinements, highest memory = 7731 / 70812

Try Ralpha Nqual Sigma-1 M(abs) CFOM Seminvariants
597829. 0.045 -0.827 0.844 1.078 0.045* +-+-- +-
Freq: 0 0 8 1 1 2 1 2 3 4 1 3 5 8 10 9 7 13 6 6 5 4 5 7 4 1 3 3 0 0 / 128
1807841. 0.043 -0.827 0.844 1.071 0.043* +-+-- +-
Freq: 0 0 15 1 3 3 3 2 4 4 3 8 14 18 25 13 14 21 13 14 9 13 8 11 10 / 256


256. Phase sets refined - best is code 1807841. with CFOM = 0.0432

Fourier and peaksearch
RE = 0.178 for 23 atoms and 487 E-values
Fourier and peaksearch
RE = 0.164 for 23 atoms and 487 E-values
Fourier and peaksearch
+*****+
+ sucrose finished at 09:27:15 Total elapsed time: 6.0 secs +
+*****+
Press any key to continue . . .
```

notes :

- The Rint and Rmerge is low (structure is solvable!)
- CFOM is low (below 0.1 structure is solved!)
- RE is low (below 20%)


Point to XSHELL



Here is the solution with a few false peaks (Q29, Q27, Q26 etc)

Move the cursor over a false peak (say Q29) and type K. Repeat for each false peak until basic structure (or fragment) is seen. Not all of the atoms may be present in this run. For example Q14 and Q26 in the example above may be misplaced. Let us delete both of them and see if one of them returns in the next run.


When finished it should look like this.



Now use the select command to choose atoms. Let us choose carbon atoms first. Move cursor to Q11 and type S: Q11 will turn blue. Repeat for Q16, Q12, Q19, Q22, Q18, Q15, Q13, Q17, Q20, Q21 and Q23 (in that order).

Now point to LABELS and GROUP LABEL... You will see


Select type C and initial number as 1 and OK. This will label the carbons (starting at Q11)



Repeat for Oxygens starting at Q5 (Q5, Q9,Q8,Q19,Q2,Q3,Q4,Q6,Q1,Q7)


Point to LABELS and GROUP LABEL... and change type to O and type OK.

Now sort the atoms. Point to ATOMS and SORT... you will see



Move all atoms to sort bin. Sort alpha numeric and insert all after <Top of List> (select <Top of List> to proceed). Type OK.

Point to Refine



point to OK and proceed.

sucrose.ins

File Edit Refine Select Atoms View Preferences Labels Render


Refine

Read instructions and data

Data: 3297 unique, 0 suppressed R(int) = 0.0535 R(sigma) = 0.0826
 Systematic absence violations: 0 Bad equivalents: 0
 wR2 = 0.5345 before cycle 1 for 3297 data and 89 / 89 parameters
 GooF = S = 3.815; Restrained GooF = 3.814 for 1 restraints
 Mean shift/esd = 2.043 Maximum = -20.204 for OSF at 10:50:19
 Max. shift = 0.085 A for C10 Max. dU = -0.016 for C7
 wR2 = 0.4858 before cycle 2 for 3297 data and 89 / 89 parameters
 GooF = S = 2.926; Restrained GooF = 2.925 for 1 restraints
 Mean shift/esd = 1.140 Maximum = -5.305 for OSF at 10:50:20
 Max. shift = 0.034 A for C12 Max. dU = -0.008 for C10
 wR2 = 0.4791 before cycle 3 for 3297 data and 89 / 89 parameters
 GooF = S = 2.837; Restrained GooF = 2.837 for 1 restraints
 Mean shift/esd = 0.382 Maximum = 1.771 for z O8 at 10:50:20
 Max. shift = 0.017 A for O8 Max. dU = 0.002 for O8
 wR2 = 0.4785 before cycle 4 for 3297 data and 89 / 89 parameters
 GooF = S = 2.815; Restrained GooF = 2.815 for 1 restraints
 Mean shift/esd = 0.196 Maximum = -0.833 for z O4 at 10:50:20
 Max. shift = 0.008 A for O4 Max. dU = 0.001 for O1
 wR2 = 0.4784 before cycle 5 for 3297 data and 2 / 89 parameters
 GooF = S = 2.810; Restrained GooF = 2.809 for 1 restraints
 R1 = 0.1924 for 2009 Fo > 4sig(Fo) and 0.2433 for all 3297 data
 wR2 = 0.4784, GooF = S = 2.810, Restrained GooF = 2.809 for all data
 R1 = 0.2447 for 1803 unique reflections after merging for Fourier
 Highest peak 5.32 at 0.4606 0.7874 0.6778 [1.43 A from C12]
 Deepest hole -0.80 at 0.1026 0.2610 0.5364 [1.32 A from O3]

OK Print

Notice the R1 is large (.1924) we are still missing an atom (the peak that is 1.43 A from C12). Point to OK



Select Q1 by moving the cursor over it and typing S. Use Group Label to Label it O11. Kill the remaining Q peaks by pointing to EDIT this Kill all Q peaks. Point to atoms and sort to resort atom list.

Point to refine again and repeat refinement.


The R1 should now be $R1 = 0.0986$. All non hydrogen atoms are now located. Find the highest Q peak by moving the cursor over Q1 (find it near C11) and reading the Peak height value in the upper left hand corner. If this is less than 2 than it is most likely a Hydrogen or residual peak. If greater than 2 then you have missed a non hydrogen atom.

Kill and Q-peaks (point to edit/Kill Q-peaks) and add hydrogen atoms.

Point to ATOMS/Hybridize ALL and then ATOMS/Calculate Hydrogens. You will see warning about HFIX 83 type OK.

Now point to refine again and repeat refinement. The $R1 = 0.0775$ now.

The final step will be the go anisotropic on all non-hydrogen atoms.



Press OK. The R1 = 0.0613


```

GooF = S = 0.873;  Restrained GooF = 0.873 for 1 restraints
Mean shift/esd = 0.018  Maximum = -0.155 for tors H8A  at 18:43:58
Max. shift = 0.009 A for H8A  Max. dU = 0.000 for C5
wR2 = 0.1897 before cycle 5 for 3297 data and 2 / 216 parameters
GooF = S = 0.873;  Restrained GooF = 0.873 for 1 restraints
R1 = 0.0613 for 2009 Fo > 4sig(Fo) and 0.1103 for all 3297 data
wR2 = 0.1897, GooF = S = 0.873, Restrained GooF = 0.873 for all data
R1 = 0.0999 for 1803 unique reflections after merging for Fourier
Highest peak 0.36 at 0.1916 0.9207 0.7949 [ 0.13 A from O7 ]
Deepest hole -0.30 at 0.4470 0.3450 0.2406 [ 1.16 A from O11 ]


*****
+ sucrose finished at 18:43:58 Total CPU time: 4.2 secs +
*****

```

Draw the structure, kill all Q-peaks (EDIT/KILL ALL Q-PEAKS) . Point to View/Rendition/Thermal Ellipsoid




Draw the thermal ellipsoid plot.



The ellipsoids should be football shaped, not disks or saucers. Non Positive Definite atoms will appear as regular shaped round gray spheres (like those seen for hydrogens and isotropic atoms).

Run the final refinement. First edit the INS file. Point to EDIT/EDIT CURRENT FILE and add

```
SIZE .2 .2 .2
TEMP -163
HTAB
CONF
WPDB -1
```



```
sucrose.ins - Notepad
File Edit Format View Help
TITL sucrose in P2(1)
CELL 0.71073 7.7576 8.7022 10.8611 90.000 102.946 90.000
ZERR 2.00 0.0018 0.0020 0.0025 0.000 0.004 0.000
LATT -1
SYMM -X, 0.5+Y, -Z
SFAC C H O
UNIT 24 44 22
SIZE .2 .2 .2
HTAB
CONF
WPDB -1
TEMP -163
L.S. 4
BOND
FMAP 2
PLAN 20
WGHT 0.081400 0.000000
FVAR 0.330420
C1 1 0.204779 0.472180 0.640813 11.000000 0.038360 0.036290 =
0.034890 -0.000750 0.004020 0.004990
AFIX 13
H1 2 0.281103 0.463850 0.725650 11.000000 -1.200000
```

Save by closing file and keep the changes. Point to refine and choose the ACTA button and start the refinement.

```

wR2 = 0.1598 before cycle 4 for 3297 data and 216 / 216 parameters
GooF = S = 1.003; Restrained GooF = 1.003 for 1 restraints
Mean shift/esd = 0.002 Maximum = 0.011 for tors H8A at 19:04:21
Max. shift = 0.001 A for H8A Max. dU = 0.000 for C1
wR2 = 0.1598 before cycle 5 for 3297 data and 2 / 216 parameters
GooF = S = 1.003; Restrained GooF = 1.003 for 1 restraints
R1 = 0.0597 for 2009 Fo > 4sig(Fo) and 0.1096 for all 3297 data
wR2 = 0.1598, GooF = S = 1.003, Restrained GooF = 1.003 for all data
R1 = 0.0992 for 1803 unique reflections after merging for Fourier
Highest peak 0.31 at 0.1959 0.9166 0.7976 [ 0.18 A from O7 ]
Deepest hole -0.30 at 0.4471 0.3452 0.2406 [ 1.16 A from O11 ]

*****
+ sucrose finished at 19:04:22 Total CPU time: 4.3 secs +
*****

```


OK Print

The final results can be examined in detail.

- 1) The R1 = 0.0597 (for Strong reflections) should be less than 0.10
- 2) The wR2 should be about 2 times R1 for all unique reflections, in this case R1=0.0992. wR2 should be less than 2 * .099 or 0.198
- 3) The GooF (1.003) should be greater than 1.0 and less than 2.0
- 4) The Maximum = 0.011 is the max. shift and should be less than 0.2
- 5) The highest peak (unassigned) is 0.31 electrons and should be less than 1.00 and the absolute value should be about equal to the absolute value of the deepest hole.


All parameters indicate that this refinement is finished.

Press ok and redraw the structure. Again kill the q-peaks. Redraw the thermal parameters. Point to label and turn off Hydrogen labels. Now move the cursor over an atom and type P. The label now will move with the cursor. Move it to a free area and click the mouse. Repeat for all atoms.



Point to render/color/to high-quality JPEG file and name the file sucrose.jpg.

Use the anti-aliasing for best pictures



Now exit xshell and point to XCIF

```

C:\WINDOWS\System32\cmd.exe
Structure Code: sucrose
[S] Change structure code
[X] Print from SHELXTL XTEXT format file
[R] Use another CIF file to resolve ? items
[C] Set compound name for tables (currently 'sucrose')
[N] Set next table number (currently 1)
[D] Set default directory for format files
[I] Crystal/atom tables from .cif
[F] Structure factor tables from .fcf
[Q] Quit

Option [R]:
Name of reference file [sucrose.pcf]: sucrose.pcf
Select data_sucrose ? (Y or N) [Y]: Y
Name of CIF file to be modified [sucrose.cif]: sucrose.cif
Select data_sucrose ? (Y or N) [Y]: Y

```

Follow the defaults and write a new Crystallographic Information File
Exit XCIF and return to SHELXTL. Point to EDIT/EDIT .CIF and edit the CIF file.

sucrose.cif - Notepad

File Edit Format View Help

```

_cell_length_a      7.7576(18)
_cell_length_b      8.702(2)
_cell_length_c      10.861(3)
_cell_angle_alpha   90.00
_cell_angle_beta    102.946(4)
_cell_angle_gamma   90.00
_cell_volume        714.6(3)
_cell_formula_units_Z 2
_cell_measurement_temperature 110(2)
_cell_measurement_reflns_used 1529
_cell_measurement_theta_min 2.69
_cell_measurement_theta_max 19.76

_expt1_crystal_description block
_expt1_crystal_colour colorless
_expt1_crystal_size_max 0.20
_expt1_crystal_size_mid 0.20
_expt1_crystal_size_min 0.20
_expt1_crystal_density_meas 0
_expt1_crystal_density_diffn 1.591
_expt1_crystal_density_method 'not measured'
_expt1_crystal_F_000 364
_expt1_absorpt_coefficient_mu 0.143
_expt1_absorpt_correction_type multi-scan
_expt1_absorpt_correction_T_min 0.9721
_expt1_absorpt_correction_T_max 0.9721
_expt1_absorpt_process_details sadabs|

```

Add crystal information

Add multi-scan And SADABS

Finally change fine-focus to normal-focus

The CIF is now ready for publication.

At this point your structure is finished. You should have three files

Sucrose.CIF, Sucrose.FCF and Sucrose.JPG.

To write a table for WORD return to XCIF. Follow the menu for crystal/atom tables. For the filename type Sucrose.RTF. For filename extension type RTA. Choose the crystal table, the atomic coordinate table, full bond lengths/angles, anisotropic table and the hydrogen atom table. Say no (N) to all other questions.

The rich text file Sucrose.RTF can now be read into any word processing program.

Example of SUCROSE

Sucrose




Table 1. Crystal data and structure refinement for sucrose.

Identification code	sucrose	
Empirical formula	C ₁₂ H ₂₂ O ₁₁	
Formula weight	342.30	
Temperature	110(2) K	
Wavelength	0.71073 Å	
Crystal system	Monoclinic	
Space group	P2(1)	
Unit cell dimensions	a = 7.7576(18) Å	α = 90°.
	b = 8.702(2) Å	β = 102.946(4)°.
	c = 10.861(3) Å	γ = 90°.
Volume	714.6(3) Å ³	
Z	2	
Density (calculated)	1.591 Mg/m ³	
Absorption coefficient	0.143 mm ⁻¹	
F(000)	364	
Crystal size	0.20 x 0.20 x 0.20 mm ³	
Theta range for data collection	1.92 to 28.37°.	
Index ranges	-10 ≤ h ≤ 10, -11 ≤ k ≤ 11, -14 ≤ l ≤ 14	
Reflections collected	8208	
Independent reflections	3297 [R(int) = 0.0535]	
Completeness to theta = 28.37°	94.3 %	
Absorption correction	None	
Max. and min. transmission	0.9721 and 0.9721	
Refinement method	Full-matrix least-squares on F ²	
Data / restraints / parameters	3297 / 1 / 216	
Goodness-of-fit on F ²	1.029	
Final R indices [I > 2σ(I)]	R1 = 0.0594, wR2 = 0.1278	
R indices (all data)	R1 = 0.1095, wR2 = 0.1543	
Absolute structure parameter	0.8(16)	
Largest diff. peak and hole	0.303 and -0.296 e.Å ⁻³	

Table 2. Atomic coordinates ($\times 10^4$) and equivalent isotropic displacement parameters ($\text{\AA}^2 \times 10^3$) for sucrose. $U(\text{eq})$ is defined as one third of the trace of the orthogonalized U^{ij} tensor.

	x	y	z	$U(\text{eq})$
C(1)	2046(5)	4721(5)	6407(4)	37(1)
C(2)	569(5)	3548(5)	6266(4)	39(1)
C(3)	-632(5)	3899(5)	7141(4)	37(1)
C(4)	-1357(5)	5506(4)	6877(4)	35(1)
C(5)	152(5)	6689(5)	6993(3)	33(1)
C(6)	3160(5)	4537(5)	5438(3)	42(1)
C(7)	1303(5)	8332(4)	8754(3)	32(1)
C(8)	2868(5)	8352(4)	9930(3)	33(1)
C(9)	4439(5)	8582(5)	9334(3)	35(1)
C(10)	3712(5)	9643(5)	8233(4)	36(1)
C(11)	-427(5)	8942(5)	8975(4)	39(1)
C(12)	4528(5)	9434(5)	7115(4)	44(1)
O(1)	1315(3)	6272(3)	6227(2)	37(1)
O(2)	2138(4)	4789(4)	4190(3)	50(1)
O(3)	1445(4)	2120(4)	6519(3)	57(1)
O(4)	-2027(4)	2773(4)	6914(3)	52(1)
O(5)	-2479(4)	5893(4)	7705(3)	42(1)
O(6)	1088(3)	6783(3)	8289(2)	32(1)
O(7)	1837(3)	9310(3)	7875(2)	35(1)
O(8)	-1210(4)	7903(4)	9701(2)	44(1)
O(9)	2958(4)	7074(3)	10737(2)	40(1)
O(10)	5891(3)	9280(4)	10216(2)	43(1)
O(11)	4601(4)	7885(4)	6732(3)	49(1)

Table 3. Bond lengths [\AA] and angles [$^\circ$] for sucrose.

C(1)-O(1)	1.460(5)
C(1)-C(6)	1.512(5)
C(1)-C(2)	1.516(5)
C(2)-O(3)	1.414(5)
C(2)-C(3)	1.504(5)
C(3)-O(4)	1.440(5)
C(3)-C(4)	1.511(6)
C(4)-O(5)	1.425(5)
C(4)-C(5)	1.543(5)
C(5)-O(1)	1.405(4)
C(5)-O(6)	1.434(4)
C(6)-O(2)	1.426(4)
C(7)-O(7)	1.409(4)
C(7)-O(6)	1.436(4)
C(7)-C(11)	1.511(5)
C(7)-C(8)	1.553(5)
C(8)-O(9)	1.408(4)
C(8)-C(9)	1.515(5)
C(9)-O(10)	1.440(4)
C(9)-C(10)	1.516(5)
C(10)-O(7)	1.448(4)
C(10)-C(12)	1.501(5)
C(11)-O(8)	1.422(5)
C(12)-O(11)	1.416(6)
O(1)-C(1)-C(6)	105.7(3)
O(1)-C(1)-C(2)	110.3(3)
C(6)-C(1)-C(2)	113.3(3)
O(3)-C(2)-C(3)	113.3(3)
O(3)-C(2)-C(1)	104.6(3)
C(3)-C(2)-C(1)	111.1(3)
O(4)-C(3)-C(2)	107.6(3)
O(4)-C(3)-C(4)	111.3(3)
C(2)-C(3)-C(4)	109.0(3)

O(5)-C(4)-C(3)	110.7(3)
O(5)-C(4)-C(5)	110.5(3)
C(3)-C(4)-C(5)	111.0(3)
O(1)-C(5)-O(6)	110.2(3)
O(1)-C(5)-C(4)	110.8(3)
O(6)-C(5)-C(4)	109.1(3)
O(2)-C(6)-C(1)	111.2(3)
O(7)-C(7)-O(6)	110.9(3)
O(7)-C(7)-C(11)	107.6(3)
O(6)-C(7)-C(11)	110.4(3)
O(7)-C(7)-C(8)	104.9(3)
O(6)-C(7)-C(8)	107.9(3)
C(11)-C(7)-C(8)	115.1(3)
O(9)-C(8)-C(9)	116.0(3)
O(9)-C(8)-C(7)	115.1(3)
C(9)-C(8)-C(7)	101.9(3)
O(10)-C(9)-C(8)	110.9(3)
O(10)-C(9)-C(10)	111.4(3)
C(8)-C(9)-C(10)	103.0(3)
O(7)-C(10)-C(12)	109.8(3)
O(7)-C(10)-C(9)	105.3(3)
C(12)-C(10)-C(9)	114.5(3)
O(8)-C(11)-C(7)	111.7(3)
O(11)-C(12)-C(10)	113.9(3)
C(5)-O(1)-C(1)	116.2(3)
C(5)-O(6)-C(7)	113.0(3)
C(7)-O(7)-C(10)	111.7(3)

Symmetry transformations used to generate equivalent atoms:

Table 4. Anisotropic displacement parameters ($\text{\AA}^2 \times 10^3$) for sucrose. The anisotropic displacement factor exponent takes the form: $-2\pi^2 [h^2 a^{*2} U^{11} + \dots + 2 h k a^* b^* U^{12}]$

	U^{11}	U^{22}	U^{33}	U^{23}	U^{13}	U^{12}
C(1)	38(2)	36(2)	35(2)	-2(2)	4(2)	5(2)
C(2)	43(2)	34(2)	37(2)	5(2)	6(2)	0(2)
C(3)	44(2)	36(2)	31(2)	0(2)	7(2)	-2(2)
C(4)	35(2)	36(2)	34(2)	-2(2)	6(2)	-5(2)
C(5)	33(2)	40(2)	26(2)	1(2)	4(2)	-3(2)
C(6)	37(2)	48(3)	41(2)	-8(2)	9(2)	2(2)
C(7)	31(2)	32(2)	30(2)	0(2)	2(2)	1(2)
C(8)	32(2)	32(2)	34(2)	0(2)	3(2)	1(2)
C(9)	30(2)	37(2)	35(2)	-7(2)	4(2)	-4(2)
C(10)	29(2)	37(2)	41(2)	-2(2)	6(2)	-6(2)
C(11)	26(2)	44(2)	45(2)	-2(2)	6(2)	3(2)
C(12)	36(2)	52(3)	42(2)	8(2)	8(2)	2(2)
O(1)	41(2)	35(2)	35(1)	-1(1)	9(1)	-3(1)
O(2)	58(2)	54(2)	39(2)	-7(1)	12(2)	-8(2)
O(3)	60(2)	40(2)	74(2)	4(2)	23(2)	8(2)
O(4)	64(2)	45(2)	49(2)	-6(2)	19(2)	-18(2)
O(5)	36(2)	48(2)	44(2)	0(1)	11(1)	-1(1)
O(6)	32(1)	31(2)	29(1)	-2(1)	-1(1)	1(1)
O(7)	25(1)	40(2)	39(1)	9(1)	3(1)	-2(1)
O(8)	38(2)	52(2)	43(2)	-5(2)	12(1)	-4(1)
O(9)	41(2)	41(2)	36(2)	6(1)	3(1)	2(1)
O(10)	24(1)	56(2)	46(2)	-9(1)	0(1)	-4(1)
O(11)	40(2)	62(2)	47(2)	-11(2)	9(2)	-2(2)

Table 5. Hydrogen coordinates ($\times 10^4$) and isotropic displacement parameters ($\text{\AA}^2 \times 10^{-3}$) for sucrose.

	x	y	z	U(eq)
H(1)	2823	4636	7273	44
H(2)	-135	3557	5374	46
H(3)	48	3826	8037	44
H(4)	-2079	5537	5990	42
H(5)	-369	7718	6718	40
H(6A)	4153	5278	5621	50
H(6B)	3665	3488	5498	50
H(8)	2748	9294	10431	40
H(9)	4812	7581	9024	42
H(10)	3879	10733	8528	43
H(11A)	-215	9940	9422	46
H(11B)	-1254	9121	8151	46
H(12A)	3839	10036	6398	52
H(12B)	5743	9856	7325	52
H(2A)	2312	5687	3961	75
H(3A)	694	1410	6435	85
H(4A)	-2358	2606	7587	77
H(5A)	-3153	6616	7391	63
H(8A)	-1537	7107	9275	65
H(9A)	3082	6269	10338	60
H(10A)	6845	8874	10142	65
H(11)	3569	7534	6508	74