

TWINNED/SPLIT CRYSTALS

CELL_NOW SAINT TWINABS

Instructor J. Reibenspies, Ph. D.
Nattamai Bhuvanesh, Ph.D.

Version 1.0.1

Start a command prompt (point START/RUN  and type cmd). Navigate to your data directory and type cell_now

CELL_NOW - Find unit-cell and matrices for twins and other problem crystals

240 reflections read from file: dpcl.p4p

read in a good P4P file. The file should have a minimum of 100 reflections and no more than 500 reflections.

Searching for vectors with $5.00 < d < 40.00$, superlattice threshold = 10.0%

Length, fom, #hits, %(0.1), %(0.2), %(0.3), components. Followed by angle/correlation to previous vectors:

....

Lots of vectors here. Ignore them unless you're insane!

The following cells would appear to be plausible, but should be checked using XPREP because they are not necessarily the conventional cells.

FOM, % within 0.2, a..gamma, volume and lattice type for potential unit-cells:

1	1.000	72.1	7.881	8.293	11.030	112.34	107.23	115.58	501.6	P
2	0.851	72.1	7.881	8.631	11.496	105.55	113.67	119.81	502.4	P
3	0.779	99.2	16.143	7.881	28.021	90.29	97.63	89.35	3533.0	P
4	0.686	99.2	16.143	7.881	30.424	89.92	114.10	90.65	3532.9	P
5	0.634	63.8	7.881	8.701	16.561	94.59	114.74	98.11	1008.9	P
6	0.607	71.2	7.881	12.493	16.143	75.91	89.35	79.65	1515.6	P
7	0.529	71.7	7.881	8.293	26.141	94.37	98.22	115.58	1507.0	P
8	0.520	71.2	7.881	8.293	26.200	78.12	80.91	64.42	1506.8	P
9	0.518	71.2	7.881	13.525	16.143	102.63	90.65	114.63	1516.0	P
10	0.504	99.2	7.881	16.143	29.144	82.48	74.02	89.35	3532.7	P
11	0.502	98.8	7.881	16.143	29.072	97.16	105.43	90.65	3533.5	P

Ignore solutions with FOM less than .5

```

1 1.000 72.1 7.881 8.293 11.030 112.34 107.23 115.58 501.6 P
2 0.851 72.1 7.881 8.631 11.496 105.55 113.67 119.81 502.4 P
3 0.779 99.2 16.143 7.881 28.021 90.29 97.63 89.35 3533.0 P

```

These three solutions stand out. #1 is the smallest cell, but only 72% of the data fits this cell. Number 3 is a much bigger cell (about 7 times) but almost all of the data fits. Solution # 2 is a perturbation of solution #1. Since I know that the Volume of the unit cell should be about 500 I am going to choose this cell. Solution # 3 is a supercell of number 1. Therefore I have chosen cell #1 with a FOM of 1.0.

Cell for domain 1: 7.881 8.293 11.030 112.34 107.23 115.58


Figure of merit: 0.593 % (0.1): 55.4 % (0.2): 72.1 % (0.3): 78.8

Orientation matrix: -0.11126408 -0.16391023 -0.08167753
0.09265278 -0.02068751 0.06747562
-0.08652867 -0.01022068 0.05087350

Percentages of reflections in this domain not consistent with lattice types:
A: 55.2, B: 54.6, C: 46.6, I: 46.6, F: 78.2, O: 67.2 and R: 67.2%

Percentages of reflections in this domain that do not have:
h=2n: 59.2, k=2n: 59.8, l=2n: 43.7, h=3n: 75.9, k=3n: 66.7, l=3n: 65.5%

RLATT color-coding employed in file: dpcn1.p4p White: indexed for first domain
Red: not yet indexed



RLATT is a program supplied by BRUKER to view the Reciprocal Lattice. The white reflections are indexed and the red are not. Notice that the red reflections are between the white layers, but NOT symmetrically placed (i.e. not 1/2 of the way between the white reflections). The distance between the white lines is 7.8Å which is the same distance between the red lines.

174 reflections within 0.250 of an integer index assigned to domain 1,
174 of them exclusively; 66 reflections not yet assigned to a domain

This line tells you how many of the original 240 reflections were indexed. In this case 174 fit the given cell and 66 do not.

In the program I have written the first cell to dpcn1.p4p

```
Cell for domain 1: 7.881 8.293 11.030 112.34 107.23 115.58
Figure of merit: 0.593 %(<0.1>): 55.4 %(<0.2>): 72.1 %(<0.3>): 78.8
Orientation matrix: -0.11126408 -0.16391023 -0.08167753
                   0.09265278 -0.02068751 0.06747562
                   -0.08652867 -0.01022068 0.05087350
Maximum deviation from integer index [0.25]:
Percentages of reflections in this domain not consistent with lattice types:
A: 55.2, B: 54.6, C: 46.6, I: 46.6, F: 78.2, O: 67.2 and R: 67.2%
Percentages of reflections in this domain that do not have:
h=2n: 59.2, k=2n: 59.8, l=2n: 43.7, h=3n: 75.9, k=3n: 66.7, l=3n: 65.5%
New cell from list <number>, reorientate <R>, accept <A> or quit <Q> [A]:
.p4p or .spin file to write domain to: dpcn1.p4p
RLATT color-coding employed in file: dpcn1.p4p
White: indexed for first domain
Red: not yet indexed
```

Next I Search for another cell that has 7.88 8.293 11.030 112.34 107.23 and 115.58 cell dimensions. I find it.

```
-----
Cell for domain 2: 7.881 8.293 11.030 112.34 107.23 115.58
Figure of merit: 0.886 %(0.1): 100.0 %(0.2): 100.0 %(0.3): 100.0
Orientation matrix: 0.10375278 0.16385944 0.08160987
                   0.05416423 0.02197985 -0.06778725
                   -0.12146208 0.00812149 -0.05056686
Rotated from first domain by 179.9 degrees about
reciprocal axis 1.000 -0.444 -0.425 and real axis 1.000 0.010 -0.002
Twin law to convert hkl from first to 1.007 0.019 -0.003
this domain (SHELXL TWIN matrix): -0.893 -1.009 0.001
                                   -0.851 -0.007 -0.998
```

The two cells are related by a rotation of 180° around the a axis (7.88Å axis)

RLATT color-coding employed in file: dpcn2.p4p
 White: indexed for first domain
 Green: current domain (but not in a previous domain)
 Red: not yet indexed

164 reflections within 0.250 of an integer index assigned to domain 2,
 66 of them exclusively; 0 reflections not yet assigned to a domain

Of the 66 remaining reflections I can fit all 66 to this new cell!

I now create a file called dpcn2.p4p which I will use in the SAINT data reduction program.

=====

Content of dpcn2.p2p

```

FILEID FRAMBO-WNT V 4.00 08/01/05 08:53:09 dpcl
SITEID Texas A&M University Bruker Analytical X-Ray Sys
TITLE dpcl
CHEM ?
CELL 7.8811 8.2935 11.0304 112.3441 107.2251 115.5824 501.623
CELLSD 0.0016 0.0017 0.0022 0.0300 0.0300 0.0300 0.251
ORT1 -0.111264080 -0.163910225 -0.081677526
ORT2 0.092652783 -0.020687509 0.067475617
ORT3 -0.086528674 -0.010220677 0.050873496
ZEROS 0.000000 0.000000 0.000000 0.0000 0.0000 0.0000
ADCOR -3.0386 -1.1271 0.0630 0.0000 -0.1217 0.3552

```

Cell number 1 with orientation matrix and zeros and corrections

```

CELL2 7.8811 8.2935 11.0304 112.3441 107.2251 115.5824 501.623
CELLSD2 0.0016 0.0017 0.0022 0.0300 0.0300 0.0300 0.251
ORT12 0.103752784 0.163859442 0.081609868
ORT22 0.054164235 0.021979848 -0.067787252
ORT32 -0.121462077 0.008121490 -0.050566856
ZEROS2 0.000000 0.000000 0.000000 0.0000 0.0000 0.0000
ADCOR2 -3.0386 -1.1271 0.0630 0.0000 -0.1217 0.3552

```

Cell number 2 with orientation matrix and zeros and corrections

```


SOURCE CU 1.54184 1.54056 1.54439 2.00000 40.00 40.00
LIMITS 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00 0.00
MORPH ?
DNSMET ?
CCOLOR ?
CSIZE ? ? ? ? ?
ADPAR 255.6000 255.7000 5.0000 512
REF05 ACHS 0 0 4 325.000 292.850 0.000 54.736 293.25 309.65 214861
189.4 0.326837 -0.269797 -0.204985

```

... etc ...


SAINT will read this file and know automatically to refine BOTH orientations. SAINT will create data files named *.mul, which we will input to TWINABS.

Start SAINT and initialize program. Now BROWSE the matrix file name and find



dpcn2.p4p. Check the cell by pointing to Cell button. If ok then delete any names in the remaining Matrix (.p4p) Filename boxes. Be sure to check the more integration files button to check ALL of the raw data sets are covered. Point to the Integrate button. Make sure that the box size is reasonable and that the

enable orientation updating is unchecked. Now point to Integrate+Sort+Global and start the integration. You will see a new set of lines in the output, the ALL, 1 and 2. These have the same meaning as those in the normal SAINT output only formatted to show both "twin" components.



```

Profile X,Y,Z spacing (deg):      0.172  0.172  0.500
Profile convolver halfwidth:     1.0    1.0    3.0
Profile convolver expansion:     1.005  1.005  2.500
Rescan threshold is top of A/D conversion range
Background pixels updated = 81.14%
Comp # File #Ref ErrX ErrY ErrZ RmsX RmsY RmsZ Inorm #Sig %<2s <Cor> %Ful XSiz YSiz ZSiz Beam
ALL 0 .001 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 0 1.20 1.20 0.60 1.000
 1 0 .001 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 0 1.20 1.20 0.60 1.000
 2 0 .001 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 0 1.20 1.20 0.60 1.000

ALL 1 .002 1 -3.07 -2.32 -0.21 3.07 2.32 0.21 991.44 5 0 0.11 40 1.20 1.20 0.60 1.000
 1 1 .002 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 1 1.20 1.20 0.60 1.000
 2 1 .002 1 -3.07 -2.32 -0.21 3.07 2.32 0.21 991.44 5 0 0.11 20 1.20 1.20 0.60 1.000

ALL 2 .003 1 0.11 -0.08 0.04 0.11 0.08 0.04 21.738 1 100 0.02 9 1.20 1.20 0.60 1.000
 1 2 .003 1 0.11 -0.08 0.04 0.11 0.08 0.04 21.738 1 100 0.02 27 1.20 1.20 0.60 1.000
 2 2 .003 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 0 1.20 1.20 0.60 1.000

ALL 3 .004 3 -2.29 0.64 0.10 2.61 1.20 0.38 1760.6 5 33 0.11 13 1.20 1.20 0.60 1.000
 1 3 .004 2 -1.48 1.24 0.32 1.60 1.41 0.39 216.09 2 50 0.10 17 1.20 1.20 0.60 1.000
 2 3 .004 1 -3.92 -0.57 -0.34 3.92 0.57 0.34 4849.6 13 0 0.11 10 1.20 1.20 0.60 1.000

ALL 4 .005 1 -1.23 -1.82 -0.38 1.23 1.82 0.38 408.82 3 0 0.26 20 1.20 1.20 0.60 1.000
 1 4 .005 0 0.00 0.00 0.00 0.00 0.00 0.00 0 0 0 0.00 0 1.20 1.20 0.60 1.000
 2 4 .005 1 -1.23 -1.82 -0.38 1.23 1.82 0.38 408.82 3 0 0.26 20 1.20 1.20 0.60 1.000

ALL 5 .006 2 -2.59 0.66 -0.05 2.60 0.93 0.06 32427 48 0 0.39 7 1.20 1.20 0.60 1.000
 1 5 .006 1 -2.37 1.31 -0.01 2.37 1.31 0.01 30041 46 0 0.02 7 1.20 1.20 0.60 1.000
 2 5 .006 1 -2.82 0.01 -0.08 2.82 0.01 0.08 34812 50 0 0.54 7 1.20 1.20 0.60 1.000

```

in the data0m._ls listing you will see that 6546 *reflections* belonged to only 1 component and 835 *reflections* belonged to both components (i.e. overlapped). In this case 9 files

```
Number of spots with 1 component(s) =      6546
Number of spots with 2 component(s) =      835 (excluded from histograms)
Number of spots with 3 component(s) =       0 (excluded from histograms)
Number of spots with 4 component(s) =       0 (excluded from histograms)
```

ending in *.mul were created. These files will be read by the next program TWINABS to generate the *hkl* file(s) we will use in the twin deconvolution.

TWINABS.

Start a command prompt (point START/RUN  and type cmd). Navigate to your data directory and type TWINABS.

```
Directory of C:\Documents and Settings\jhr6675\Desktop\class\twin\work
05/22/2006 10:49a      164,223 data01.mul
05/22/2006 10:52a      163,709 data02.mul
05/22/2006 10:54a      160,368 data03.mul
05/22/2006 10:57a      278,331 data04.mul
05/22/2006 11:00a      255,972 data05.mul
05/22/2006 11:03a      264,967 data06.mul
05/22/2006 11:06a      267,280 data07.mul
05/22/2006 11:10a      284,242 data08.mul
05/22/2006 11:14a      272,420 data09.mul
05/22/2006 11:15a      2,111,512 data0n.mul
                4,223,024 bytes free
    10 File(s)          3,300,950,016 bytes free
    0 Dir(s)

C:\Documents and Settings\jhr6675\Desktop\class\twin\work>twinabs

TWINABS - Bruker Nonius scaling and corrections for twinned crystals - 01.05

Maximum number of reflections allowed <500000>:
Enter listing filename [twin.abs]:

Laue group numbers:

[1] -1                      [8] -3m (rhombohedral axes)
[2] 2/m (Y unique)         [9] -3in (Z unique)
[3] mm                      [10] -3m1 (Z unique)
[4] 4/n (Z unique)         [11] 6/n (Z unique)
[5] 4/mmm (Z unique)       [12] 6/mmm (Z unique)
[6] -3 (rhombohedral axes) [13] m3
[7] -3 (Z unique)          [14] m3m

[0] to write list of equivalent indices for Laue/point groups to listing file
Enter Laue group number [2]: 1

Treat Friedel opposites as equivalent for parameter refinement (Y or N) ?
Answering "N" halves the data to parameter ratio and is not recommended
unless you have a high redundancy and know what you are doing [Y]:

Read reflection files written by EVALCCD with extension .sam specified) or by
SAINT (extension .mul, default if no extension). Either individual files for
each scan or a single merged file may be read. It is important that all files
are from the same crystal and that reflections have been indexed consistently,
i.e. that the orientation matrices are similar (no rows with signs reversed)!
```

Choose the correct Laue group and answer Y to the next question. Enter the file name for the first data set, in this case it is data01.mul. Now enter the rest one by one and end with a /.

It is better to use the unmerged data sets at this point. You will

need the unmerged data later. Fit the data to a single component (in this case component 2) In this case use the defaults and go to the refinement step

```
Enter filename </ if no more> [ 1]: data01
Enter filename </ if no more> [data02.mul]:
Enter filename </ if no more> [data03.mul]:
Enter filename </ if no more> [data04.mul]:
Enter filename </ if no more> [data05.mul]:
Enter filename </ if no more> [data06.mul]:
Enter filename </ if no more> [data07.mul]:
Enter filename </ if no more> [data08.mul]:
Enter filename </ if no more> [data09.mul]:
Enter filename </ if no more> [data010.mul]: /

    2 twin components present

Mean and maximum errors in direction cosine check function =  0.001  0.005
The mean error should not exceed 0.008, and is usually caused by matrix
changes during data processing.

Maximum 2-theta = 124.89 deg.   Approximate wavelength = 1.54178 Angstroms
=====
PART 1 - Refinement of parameters to model systematic errors

3273 data ( 1201 unique ) involve component 1 only, mean I/sigma 22.1
3273 data ( 1205 unique ) involve component 2 only, mean I/sigma 21.2
 835 data ( 400 unique ) involve 2 components, mean I/sigma 27.8

Enter N to fit component N and apply to all data, e.g. when one component
dominates; -N for one parameter set to fit equivalent singles and composites
that contain at least one contribution from domains 1..N; or 0 to fit each
component separately [-2]: -2
```

The refinement will model the absorption surface. The refinement went well so accept the results. Go to the outlier rejection step

```

Refinement of a single parameter set to fit all single and composite
reflections that contain at least one contribution from domains 1.. 2
-----
  4711 single and   550 overlapped reflections used for parameter fitting
Effective data to parameter ratio =   0.96

R<int> =  0.0585 (selected reflections only, before parameter refinement)

Cycle  R<incid>  R<diffr>  Mean wt.
  1     0.0536   0.0512   0.9379
  2     0.0504   0.0495   0.9414
  3     0.0497   0.0495   0.9420
  4     0.0497   0.0495   0.9422
  5     0.0496   0.0495   0.9423
  6     0.0497   0.0495   0.9422
  7     0.0497   0.0496   0.9422
  8     0.0497   0.0496   0.9422
  9     0.0497   0.0496   0.9421
 10     0.0498   0.0497   0.9421
 11     0.0498   0.0497   0.9421
 12     0.0498   0.0497   0.9420
 13     0.0499   0.0498   0.9420
 14     0.0499   0.0498   0.9419
 15     0.0499   0.0498   0.9419

R<int> =  0.0498 (selected reflections only, after parameter refinement)
Repeat parameter refinement <R> or accept <A> [A]: a

```

and use the defaults. The g is predicted to be 0.085 which is good and should be accepted. Do not accept anything over 0.2.

Statistics for singles for twin component 1									
Run	2-theta	R<int>	Incid. factors	Diffr. factors	K	Total I>2sig(I)			
1	-35.0	0.0504	2.805	-3.063	0.897	-1.063	0.752	249	221
2	-35.0	0.0473	2.463	-2.856	0.868	-1.038	0.751	262	234
3	-35.0	0.0399	2.903	-3.234	0.880	-1.113	0.765	273	244
4	-90.0	0.0359	0.366	-0.417	0.897	-1.202	0.604	391	334
5	-90.0	0.0316	0.351	-0.416	0.863	-1.042	0.613	380	324
6	-90.0	0.0353	0.311	-0.383	0.880	-1.111	0.580	385	322
7	-90.0	0.0267	0.296	-0.358	0.863	-1.054	0.518	368	316
8	-60.0	0.0407	0.875	-1.058	0.864	-1.427	0.627	449	409
9	-60.0	0.0452	0.802	-0.947	0.879	-1.002	0.622	445	409

Statistics for singles for twin component 2									
Run	2-theta	R<int>	Incid. factors	Diffr. factors	K	Total I>2sig(I)			
1	-35.0	0.0466	2.798	-3.053	0.898	-1.062	0.856	245	209
2	-35.0	0.0399	2.463	-2.855	0.867	-1.043	0.766	264	231
3	-35.0	0.0475	2.901	-3.243	0.880	-1.184	0.784	274	242
4	-90.0	0.0305	0.365	-0.417	0.898	-1.326	0.581	402	342
5	-90.0	0.0294	0.351	-0.416	0.865	-1.056	0.583	371	312
6	-90.0	0.0374	0.311	-0.382	0.880	-1.104	0.615	382	314
7	-90.0	0.0335	0.298	-0.358	0.865	-1.059	0.572	380	307
8	-60.0	0.0471	0.878	-1.056	0.862	-1.417	0.654	438	389
9	-60.0	0.0424	0.792	-0.947	0.882	-1.008	0.691	443	393

Statistics for all composite reflections									
Run	2-theta	R<int>	Incid. factors	Diffr. factors	K	Total I>2sig(I)			
1	-35.0	0.0368	2.800	-3.063	0.907	-0.968	0.856	63	52
2	-35.0	0.0234	2.540	-2.846	0.888	-1.059	0.559	50	47
3	-35.0	0.0201	3.133	-3.245	0.880	-1.058	0.485	31	29
4	-90.0	0.0331	0.367	-0.414	0.907	-1.016	0.554	132	118
5	-90.0	0.0524	0.351	-0.416	0.880	-1.059	0.645	118	106
6	-90.0	0.0313	0.317	-0.382	0.884	-1.036	0.550	120	108
7	-90.0	0.0420	0.299	-0.358	0.863	-1.060	0.651	139	120
8	-60.0	0.0425	0.900	-1.058	0.863	-1.207	0.730	92	87
9	-60.0	0.0371	0.794	-0.947	0.909	-0.970	0.637	77	74

Statistics for all single and composite reflections									
Run	2-theta	R<int>	Incid. factors	Diffr. factors	K	Total I>2sig(I)			
1	-35.0	0.0463	2.798	-3.063	0.897	-1.063	0.811	557	484
2	-35.0	0.0400	2.463	-2.856	0.867	-1.059	0.745	576	512
3	-35.0	0.0403	2.901	-3.245	0.880	-1.184	0.761	578	514
4	-90.0	0.0331	0.365	-0.417	0.897	-1.326	0.588	925	795
5	-90.0	0.0347	0.351	-0.416	0.863	-1.059	0.604	869	744
6	-90.0	0.0356	0.311	-0.383	0.880	-1.111	0.592	887	744
7	-90.0	0.0324	0.296	-0.358	0.863	-1.060	0.558	887	746
8	-60.0	0.0436	0.875	-1.058	0.862	-1.427	0.650	979	885
9	-60.0	0.0427	0.792	-0.947	0.879	-1.008	0.656	965	876

su = K * Sqrt[sigma^2<I> + <g<I>>^2] where sigma<I> is estimated by SAINT

The statistics for the data set are good. Note the low R(int) and the similar K values. This indicates an acceptable correction. Keep it and continue.

Output reflections, here we choose to output a standard *hkl* file with only one component, which we name DPA_HKLF4.hkl and a non-standard *hkl* file with both twin components (1 and 2) as

```
PART 3 - Output Postscript diagnostics and corrected data

Write Postscript diagnostic file <Y or N> [Y]: n
It is now possible to generate a HKLF 4 format file with crude intensities
for structure solution or a more accurate HKLF 5 format file for refinement.
Repeat <R>, write .hkl file <4 or 5>, or quit <Q> [4]: 4
Average equivalent reflections <Y or N> [Y]: Y
Enter name of output .hkl file [twin.hkl]: DPA_HKLF4.HKL
Mu*r of equivalent sphere for additional spherical absorption correction.
Enter <CR> if none:
Average Friedel opposites in output file <Y or N> [Y]: Y
Use component N only <recommended>; 0 to use all [1]: 1

    1471 Corrected reflections written to file DPA_HKLF4.HKL

Reflections merged according to point-group -1
Ratio of minimum to maximum apparent transmission: 0.604371

It is now possible to generate a HKLF 4 format file with crude intensities
for structure solution or a more accurate HKLF 5 format file for refinement.
Repeat <R>, write .hkl file <4 or 5>, or quit <Q> [5]: 5
Average equivalent reflections <Y or N> [Y]: Y
Enter name of output .hkl file [twin.hkl]: DPA_HKLF5.HKL
Mu*r of equivalent sphere for additional spherical absorption correction.
Enter <CR> if none:
Average Friedel opposites in output file <Y or N> [Y]: Y
Output only reflections containing component N <0 to output all> [1]: 0

    3193 Corrected reflections written to file DPA_HKLF5.HKL

Reflections merged according to point-group -1
Ratio of minimum to maximum apparent transmission: 0.603891

It is now possible to generate a HKLF 4 format file with crude intensities
for structure solution or a more accurate HKLF 5 format file for refinement.
Repeat <R>, write .hkl file <4 or 5>, or quit <Q> [Q]:
```

DPA_HKL5.hkl. Notice for the first file I asked for component 1 when prompted and for the next file I asked for both. Also the first file is in HKLF 4 format and the 2nd file is in HKLF 5 format. We will use them both later.

You now have two *hkl* files. Use the HKLF4 file to solve the structure and refine it. If the solution and refinement look good that is all you may want to do.

I solved and refined the structure using the HKLF 4 data set in the normal way. The residuals refined to

_refine_ls_R_factor_all	0.0561
_refine_ls_R_factor_gt	0.0488
_refine_ls_wR_factor_ref	0.1832
_refine_ls_wR_factor_gt	0.1711
_refine_ls_goodness_of_fit_ref	1.076
_refine_ls_restrained_S_all	1.076
_diffn_reflns_theta_full	62.12
_diffn_measured_fraction_theta_full	0.899

This is perfectly acceptable, but we can do better if we refine the HKLF5 data. In this case insert the line **BASF 0.5** before the FVAR line and change the **HKLF 4 line to HKLF 5**. Now refine the file with the DPA_HKLF5.hkl data. The new residuals are

_refine_ls_R_factor_all	0.0459
_refine_ls_R_factor_gt	0.0397
_refine_ls_wR_factor_ref	0.1206
_refine_ls_wR_factor_gt	0.1130
_refine_ls_goodness_of_fit_ref	1.056
_refine_ls_restrained_S_all	1.056
_diffn_reflns_theta_full	62.55
_diffn_measured_fraction_theta_full	0.962

The R(all) decreased more than a percent while the wR(all) decreased more than 6 %. Also the percent coverage increased from ~90% to ~96%.

The BASF which is the refined component scale is ~ 0.49 or a 50/50 mixture of twin components.