

Daily Safety Check: room 2407 and rm 2409

Date : _____

- D8 Focus Bragg-Brentao
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Specimen Holder OK
 - Enclosure OK
- D8 Vario
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Specimen Holder OK
 - Enclosure OK
 - No unsafe conditions
- SAXS RAG long
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Vacuum pump OK
 - Enclosure OK
 - No unsafe conditions
- SMART 1000 Workstation 1
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Beam Stop in Place
 - CEU ok (not overheating)
 - Enclosure OK
 - No unsafe conditions
 - LT tanks OK
 - LT attachment normal
- SMART 1000 Workstation 2
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Beam Stop in Place
 - CEU ok (not overheating)
 - Enclosure OK
 - No unsafe conditions
 - LT tanks OK
 - LT attachment normal
- GADDS
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Beam Stop in Place
 - Hi-Star OK
 - Enclosure OK
 - No unsafe conditions
 - LT tanks OK
 - LT attachment normal
- APEX Workstation 3
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Beam Stop in Place
 - APEXII water OK
 - Enclosure OK
 - No unsafe conditions
 - LT tanks OK
 - LT attachment normal
- SAXS Short
 - X-rays ON and Normal
 - SHUTTER OK
 - Water Flow OK
 - Vacuum pump OK
 - Enclosure OK
 - No unsafe conditions
- General Room Safety
 - No electric extension cords
 - Lighting OK
 - SHARPS properly stored
 - Chemicals stored or labelled properly
 - Chemical Showers OK (test)
 - No unsafe conditions
 - No food or drink in lab area
 - Evacuation routes clear
 - Fire extinguishers clear
 - Electrical Panels clear
 - Water Traps filled
 - Warning Signs OK
 - No Loose Cords on Floor
 - Gas cylinders secured